



Analysis of Pharmacologically Relevant Compounds using GC/MSD – EI/PCI/NCI

Compendium of Applications



Agilent Technologies

**Analysis of
Pharmacologically Relevant Compounds
using
GC/MSD – EI/PCI/NCI**

Compendium of Applications

Preface

The combination of gas chromatography with mass spectrometry has been successfully used for decades. The now standard and universal technique of electron impact ionization (EI) has provided, with respect to sensitivity and compound identification, extraordinary performance and the volume of spectral reference data has continued to increase over the years. Commercially available mass spectral libraries now contain more than 350,000 entries. Instrument sensitivity has improved such that sufficient sensitivity is available in full scan mode for determinations at analyte concentrations in the range of 10pg/μl to 1pg/μl. However the EI technique still leaves some analytical demands unfulfilled. For example, if the EI spectrum shows little or no definitive information about the analyte molecule, such as molecular weight, or, in samples where matrix interferences convolute the spectral information due to insufficient selectivity, then successful detection can decrease drastically.

For such cases, chemical ionization techniques (CI) are available which enhance and extend the use of mass spectrometry. The positive chemical ionization (PCI) mode can usually lead to molecular weight information since protonated and adduct species of sufficient intensity can provide distinctive information about the

correct molecular mass. For compounds with high electron affinity, the electron capture negative ion chemical ionization (ECNI or NCI) mode offers selectivity and the utmost sensitivity even in difficult matrices. Detection limits in the femtogram range can often be easily achieved. Stable instrumental parameters can insure high reproducibility over a wide concentration range. Importantly, the new instrumental platforms provide user-friendly CI operational features, such as autotuning, that make CI similar to EI in convenience.

With CI techniques, mass spectrometry becomes a powerful, unique, problem-solving tool. CI is the most creative MS technique. After becoming familiar with some criteria and the effects of certain instrumental parameters, the user will rapidly achieve positive results, and the usefulness and fascination of CI will become readily apparent. To assist and convey such an experience is the aim of this compendium.

*I would like to thank
Dr. Harry Prest, Senior
Application Chemist,
Agilent Technologies,
Chemical Solution Business
Division, Palo Alto, USA,
for his helpful comments
and advice in preparing this
applications compendium.*

H.-Jürgen Schulz
February, 2002

- 1. Introduction
- 1.1 Distinctions between EI & CI
- 2. Positive Chemical Ionization (PCI)
- 2.1 Ion Source Configuration
- 2.2 PCI – Reactions of different Reagent Gases
- 2.3 Examples of EI and PCI
- 3. Electron Capture Negative Ionization (ECNI)
- 3.1 Optimizing ECNI Analysis
- 4. Practical Operating Advice
- 4.1 Gas Chromatography
- 4.2 PCI & NCI Conditions
- 4.3 Derivatization
- 5. Instrumentation
- 6. Literature

1. Introduction

The technique of chemical ionization mass spectrometry (CIMS) was pioneered by Munson and Field in 1966. CIMS can be considered an alternative ionization approach to electron impact (EI) and offers the opportunity to determine the molecular weight of the analytes and in selected cases, CIMS measurements show very high selectivity and outstanding sensitivity.

1.1 Distinctions between EI & CI

In EI mode, relatively high energy electrons (70 eV) collide with analyte molecules producing positive ions and other species. The fragmentation process, executed under constant conditions, is well understood and the positive ion fragmentation pattern, which is the EI mass spectrum of the analyte, is used for compound identification. Whereas EI is a direct energy transfer process with electron kinetic energy deposited directly in an analyte molecule, CI is an indirect process involving an

intermediate chemical agent. This is particularly true in positive chemical ionization (PCI). In PCI, the ion source is filled with a reagent gas which is ionized to create reagent ions which react with the analyte. The “interesting” reaction products are positive ions which are collected and measured in PCI. Negative ions are also formed by this chemical process and this is a form of NCI. The gas filling the source can also be used to “buffer” or “thermalize” electrons. These “slow” electrons can be captured very efficiently by analyte molecules to form negative ions which is called electron capture negative ionization (ECNI). In casual usage, when the resulting ions being measured are positive, it is referred to as PCI and when negative, NCI.

2. Positive Chemical Ionization (PCI)

PCI mode is preferred in cases that the EI mass spectrum contains mostly low mass-to-charge ratio fragments and therefore little or no information about the molecular weight of the analyte. The success of the technique is strongly dependent on the choice of the applied reagent gas. The literature documents a multitude of examples for structure elucidations with PCI techniques by using different reagent gases, in addition to the elementary use of PCI for molecular weight. PCI sensitivity is comparable to EI sensitivity but can offer significant improvements in actual samples depending on the choice of CI reagent gas and nature of the interferences. In terms of the signal-to-noise of the base peaks, PCI often provides a better result than EI and therefore is well suited to selected ion monitoring (SIM) acquisitions.


2.1 Ion Source Configuration

The CI ion source resembles the EI source but is designed to have an ionization chamber (< 1ml volume) where the reactions take place that is much more enclosed. The filament generating the electrons is


positioned just outside the chamber. The electrons emitted from the filament into the ion chamber are accelerated to between 100 and 200 eV for optimal penetration of the reagent gas. The electron entrance orifice is small to keep the chamber tight and reagent gas pressure high. The reagent gas enters the ionization chamber via the GC/MSD interface. Because the amount of reagent gas in the source completely overwhelms the analyte, the reaction of the analyte is very efficient. The shape of the source is designed that the partial pressure in the outer chamber is about 10^{-5} to 10^{-6} torr in order to maintain a collision free path for the ions to reach the analyser. Mass spectrometry using such a source design is named as High Pressure Mass Spectrometry (HPMS).

2.2 PCI reactions of different Reagent Gases

The most frequently used reagent gases are methane, iso-butane and ammonia. The reagent gas is ionized by electrons entering the ionization source and, because the pressure of the reagent gas is high, a number of reactions occur. The principal methane reactions are:



This stepwise description of the gas phase reactions is for clarity only and in fact all the processes are proceeding simultaneously. Of special interest is the CH_5^+ ion which is a strong proton donor and can form a protonated molecule with an analyte. Additionally methane forms characteristic molecular adduct ions :



Such processes are described by the term “Proton Affinity” (PA). PA is defined as the thermochemical


ability of the reaction partners to transfer protons. In order to generate a protonated molecule, the PA of the analyte must be greater than that of the reagent gas ion. As the difference between the PA of the analyte and the PA of the reagent gas ion increases the amount of fragmentation of the analyte increases. Referring to literature, the PA of organic compounds is between 180kcal/mol to 240kcal/mol. In most cases, the PA of the compound of interest is unknown and successful PCI requires experimentation with reagent gases in each case.

Reagent Gas Ions and their Proton Affinities (kcal/mol):


Hydrogen	H_3^+	101
Methane	CH_5^+	132
iso-Butane	$C_4H_9^+$	196
Ammonia	NH_4^+	204

Ammonia has a high PA value which is closest of the common gases to the PA of most organic molecules so it often provides better differentiation (more selectively) between the sample matrix and the analyte. Non-polar compounds or matrices consisting mostly of C and H, are less readily ionized in comparison to more polar compounds. For example, carboxylic esters, such as the phthalates (Agilent Technologies Application Note 5988-2244EN) show a high affinity for ammonia because of the polar ester-linkages (COO groups).

Formation of Ammonia Reagent Gas Ions:




Typical ammonia reactions with an analyte (M)



Analytes suited for ammonia reaction generate the characteristic $[M+NH_4]^+$ adduct ion, some times the protonated molecule and most often both.

The choice of the applied reagent gas determines the fragmentation behavior of the analyte and consequently the result of the PCI measurement. PCI processes generating little fragmentation are referred to as “soft ionization” and the corresponding reagent gases are called “soft reagent gas”. Conversely, reactions creating a great deal of fragmentation are named “hard ionization” and the gas, a “hard reagent gas”.

In some cases the Hydride Abstraction Reaction is of interest, and often occurs with long chain alkanes or compounds containing long chain alkyl groups:



In these cases the PA of the reagent gas ion is greater than the PA of the analyte.

2.3 Examples of EI/PCI

Figure 1 shows an example of the differences in response for the different ionization processes and PCI reagent gases for methyl palmitate. In full scan mode the total ion response follows $EI > PCI/CH_4 > PCI/NH_3$. The spectra reflects the characteristics of each reaction, Figure 2










Figure 1. Total ion chromatograms of methyl palmitate at 600pg injected, in EI (top), PCI-CH₄ (middle) and PCI-NH₃ (bottom)

Figure 2. Full Scan Spectra referring to Figure 1

3. Electron Capture Negative Ionization (ECNI)

The importance of the ECNI technique becomes readily evident in view of the selectivity of the detection of suitable analytes at very low concentration levels (ppt). In addition, the typical sample matrix interferences are generally suppressed in ECNI MS mode resulting in very high signal-to-noise values. Suitable analytes for NCI, have a high electron capture capacity or high electron affinity (EA). The expression “chemical ionization” is not applicable because the important part of the reaction is achieved by capture of low energy electrons – “thermal” electrons – by the analyte so there are no ion molecular reactions involved in the initial ion creation process. Thermal electrons are generated by collision of electrons emitted from the filament with buffer gas (e.g. methane) molecules located at high pressure in the ionization chamber of the source:



Analytes with high EA form stable molecular anions, M^- , and show a simple spectrum. (See Figure 3). This process is resonant electron capture and commonly called NCI. Electrons with an energy potential of about 15 eV lead to dissociative reactions:



Mass spectra resulting from these reactions exhibit more fragmentation and less sensitivity compared with the ECNI process. Though it is hard to predict which analytes are suitable for NCI some rules of thumb are used. Good results are often obtained with compounds already successfully analyzed by GC/ECD. ECNI candidates contain multiple halogen groups, or nitro groups, or double bond and/or conjugated structures, and/or hetero atoms. Polar compounds which are suitable for derivatization can be modified with perfluoro reagents in order to integrate a group with a high EA into the




Figure 3. Mass spectra of tetrahydrocannabinol as the trifluoroacetic acid anhydride derivative (nominal mol. wt. 410 g/mol) in EI (top) and ECNI-CH4 (bottom)

molecular structure and improve chromatography.

3.1 Optimizing ECNI measurements

In ECNI mode some ion source operating parameters can be modified to improve sensitivity.

Buffer Gas

Raising the flow amount of the buffer gas, which influences the collision process, will sometimes lead to an increase in response. Such optimization refers to the term “High Pressure Electron Capture Mass Spectrometry” (HPECMS). The quality and purity of the buffer gas are of importance. The gas purity should be in the range of 99.95% (3.5) to 99.995% (4.5). In addition to methane, ammonia and carbon dioxide are frequently applied. Oxygen and water are very efficient at collecting electrons and can suppress analyte response and great care must be taken that they are excluded.

Ion Source Temperature

The ionization chamber temperature dramatically affects the reaction yield and the fragmentation behavior and

consequently, the sensitivity. Low ion source temperature favors the electron capture process. As a consequence of the operation of the source (hot filament) and the GC (hot carrier gas), the lowest practical source temperature is about 150°C. Considering the analyte elution temperature, the ion source body sometimes is a “cold spot” and might be a reason for tailing peak shape and matrix effects.

Tuning Parameters

The tuning program controls the parameters of electron energy (eV) and emission current (μA). Electron energy (EE) has some influence on the mobility of the electrons and their penetration efficiency in the ionization chamber. A high EE value is usually advantageous. Emission current (EC) is related to the amount of electrons emitting from the filament. Increasing the EC value leads to response additional but filament lifetime has to be considered. Higher ECs lower filament life and the accuracy of the filament position above the electron entrance slit may become degraded due to deformation.

4. Practical Hints

The following hints are based on the experience of the author and are to be considered as advice and do not represent a guarantee in all circumstances or applications.

4.1 Gas Chromatography

Chemical Ionization requires the same GC criteria as applied in EI mode. Carrier gas purifiers are highly recommended. Also an air-water check should be performed prior to switching to CI. Before developing a CI method, all GC operating parameters such as injection technique and the capillary column should be tested in EI mode in order to avoid any kind of sample discrimination. Thermally labile compounds are analysed using on-column or PTV injection systems which permit cryo-focussing and sample enrichment. Splitless injection executed in pulsed-pressure mode

has the advantages of almost quantitative sample transfer from the injection liner into the column and the reduced residence time of the sample in the injector inlet. The liner type has to be appropriate to the injection technique. In this compendium, for splitless injection a deactivated double-taper liner (Agilent Part Nr. 5181-3315) was used. The choice of the column is always related to the analytical problem, however, for the majority of the analytes documented here, the MSD standard column (HP-5ms, 30m x 0.25mm x 0.25µm, Agilent Part Nr. 19091S-433) was adequate. Sample matrix interferences may suggest use of another phase and /or GC oven program.

4.2 PCI and NCI Conditions

For both PCI and NCI techniques, the MSD software includes sophisticated autotuning programs. These autotune programs assist adjusting the flow of reagent gas and the other parameters necessary for successful CI operation. These autotune values are a good starting point. The value for Electron Multiplier Voltage (EM Voltage) normally needs to be increased in over autotune values by about 400 V. In PCI mode, the results are primarily dictated by choice and pressure of the reagent gas, altering autotune parameter values has almost no benefits. In ammonia PCI, increasing the pressure typically increases the formation of adduct ions. In NCI mode the previously discussed tune parameters are worth adjusting in order to improve sensitivity. In any case frequent tuning should be avoided, especially in ECNI mode where the residual tuning gas can increase the background for several hours due to the extreme sensitivity in ECNI. The analyte amount (absolute amount onto the column) recommended for method development in scan mode are 10ng for PCI and 1ng or less for NCI.

4.3 Derivatization

Polar, chromatographically difficult compounds frequently need derivatization before analysis. In NCI mode, derivatized analytes show increased sensitivity. The following recommendations refer

to the literature references and to the author's lab experiences. Derivatization reagents and reaction (incubation) criteria are matter of choice and may be varied in order to improve reaction yield and sensitivity. For a reaction vial, the "High Recovery Vial", 1.5ml volume, conical bottom autosampler vial (Agilent Part Nr. 5182-3454) is recommended.

Care must be taken to insure that all solvents used for derivatization are free of water. The solution containing the sample that is to be derivatized is evaporated (blown-down) with purified dry nitrogen introduced into the vial by means of a capillary steel tube (1/16" o.d.) or glass pipet. The tube's opening is positioned some millimeters above the liquid surface and a gentle gas flow (checked prior to placing over the reaction vial) which forms only a small depression on the surface, is applied until the solvent has completely evaporated. The derivatization reagent is added to the dry residue, the vial sealed, and allowed to react for a specific length of time, at a regulated and usually elevated temperature. Depending on the chemical nature of the derivatizing reagent (i.e., capable of degrading the capillary column phase), it may have to be removed by nitrogen blow-down as described above and reconstituted in an appropriate solvent before the sample is injected. Sometimes an adequate dilution must be prepared. According to experience, most derivatives are unstable, even if they are stored at low temperature. In most cases, when analyzing a sample instead of a standard, derivatization happens to analytes and to any reactive matrix compounds present. Such byproducts can complicate the chromatogram and spectral analysis and should be minimized or eliminated by further sample preparation steps or improved chromatography.

Chemically aggressive reagents will deteriorate the stationary phase of the column, especially when applying splitless or on-column injection. With split injection, care should be taken to avoid corrosion of the gas tubing

connected to the injector and that the the gas regulation module (EPC module) is not damaged.

5. Instrumentation

The documented applications were executed with the Agilent Technologies GC/MSD System:

- Gas Chromatograph 6890plus, split/splitless and On Column Injector, Autosampler 7673
- Mass Spectrometer MSD 5973N, CI Option
- HP Kayak XA, ChemStation Software Vers. G1701CA

As reagent gas or buffer gas methane (4.5) and ammonia (3.5 or 4.0), Linde Gas AG, were used. All gas supply tubes were of stainless-steel material and a gas purifier (Agilent Part Nr. 1999-80410 used for **only** methane) were installed between the MSD and the gas bottles. For ammonia, an appropriate, corrosion resistant gas regulator was used. The ammonia gas stainless steel supply tube was coiled and the pressure was adjusted to approximately 7 psi (0,5 bar) in order to avoid generating droplets. (For more information on operating in ammonia, refer to Agilent Technolgies Technique Brief Nr. 5968-7844E).

6. Literatur

"Chemical Ionization Mass Spectrometry", 2nd Edition, Alex G. Harrison, CRC Press, ISBN 08493-4254-6

"Introduction to Mass Spectrometry", Chapter Six, 2, J. Throck Watson, Raven Press, New York

"High Pressure Electron Capture Mass Spectrometry", W. B. Knighton, L. J. Sears and E. P. Grimsrud, Mass Spectrometry Reviews, 1996, 14, 327-343

"Handbook of Analytical Derivatization Reactions", D. R. Knapp, John Wiley & Sons, ISBN 0-471-03469-X

"Handbook of Derivatives for Chromatography", 2nd Edition, K. Blau and J. Halket, John Wiley & Sons, ISBN 0-471-92699-X

"Silylating Agents", Fluka, ISBN 3-905617-08-0

Contents

Component	Summary
A	
Acepromazine 7	
Alprazolam 19	
Amobarbital 9	Barbiturates 9
B	Amobarbital
Barbital 9	Barbital
Benzoylcegonine 25	Butethal
Bromazepam 19	Pentobarbital
Butethal 9	Secobarbital
C	
Chloramphenicol 29	Benzodiazepines 19
Chlorphenoxamine 33	Alprazolam
Chlorprothixene 35	Bromazepam
Cholesterol 95	Diazepam
Cimaterol 37	Flunitrazepam
Clenbuterol 41	Triazolam
Cocaine 47	
Codeine 49	Nitroimidazoles 79
D	Dimetridazole
Diazepam 19	Metronidazole
Dimethindene 53	Ronidazole
Dimetridazole 79	
Diphenhydramine 55	Steroides 95
E	Cholesterol
Estradiol 95	Estradiol
Estrone 95	Estrone
F	Testosterone
Flunitrazepam 19	
L	
Lidocaine 57	
M	
Mabuterol 59	
MDA 63	
(Methylenedioxyamphetamine)	
Mepivacaine 67	
Methadone 69	
Metronidazole 79	
Morphine 71	
Nalorphine 75	
O	
Orphenadrine 83	
P	
Pentobarbital 9	
Phenylbutazone 85	
Promethazine 87	
Propionylpromazine 89	
(Combelen)	
R	
Ractopamine 91	
Ronidazole 79	
S	
Secobarbital 9	
T	
Testosterone 95	
Tetrahydrocannabinol 101	
(THC)	
Tetrahydrocannabinol	
Carboxylic Acid 105	
(THCCOOH)	
Triazolam 19	
Z	
Zearalenone 109	

Acepromazine

CAS-Nr. 61-00-7

Molecular formula: C₁₉H₂₂N₂OS

GC-Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3 min) - 20°C/min to

300°C (4 min)

MS-Parameter

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1 ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1 ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 11.60min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 18/1


PCI/NH₃ Scan: > 25/1



EI Spectrum, Acepromazine: *m/z* 326; M⁺



PCI/CH₄ Spectrum, Acepromazine: *m/z* 327, 355, 367; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃ Spectrum, Acepromazine: *m/z* 327; [M+H]⁺

Barbiturates

Amobarbital CAS-Nr. 57-43-2
Molecular formula: C₁₁H₁₈N₂O₃
Barbital CAS-Nr. 57-44-3
Molecular formula: C₈H₁₂N₂O₃
Butethal CAS-Nr. 77-28-1
Molecular formula: C₁₀H₁₆N₂O₃
Pentobarbital CAS-Nr. 76-74-4
Molecular formula: C₁₁H₁₈N₂O₃
Secobarbital CAS-Nr. 76-73-3
Molecular formula: C₁₂H₁₈N₂O₃

GC-Parameter

Column: HP-5ms
Agilent Part Nr. 19091S-433
30m x 0.25mm x 0.25µm
Carrier Gas: Helium
Flow: 0.7ml/min, 30cm/sec
Mode: Constant Flow
Injection: Pulsed splitless, 250°C
Oven Temperature Program
60°C (1min) – 20°C/min to 180°C
10°C/min to 300°C

MS-Parameter

Mode: EI – SCAN
Tune: Atune
Temperatures:
Source 230°C, Quad 150°C
Mode: PCI/CH₄ – SCAN
Reagent Gas: Methane
Flow (Setting): 1ml/min (20)
Tune: PCI-Methane Autotune
Temperatures:
Source 250°C, Quad 106°C
EM Voltage: Tune + 400V
Mode: ECNI/CH₄ – SCAN/SIM
Buffer Gas: Methane
Flow (Setting): 2ml/min (40)
Tune: ECNI-Methane Tune File
Temperatures:
Source 150°C, Quad 106°C

Remarks


Derivatization with Pentafluorobenzylbromide (PFBB)
To 100µl of the Barbiturate
Standard (SIGMA D 3155),

concentration 20ng/µl each, diluted in ethylacetate, 10µl of the derivatization reagent and 10µl of triethylamine are added and the mixture is incubated for 60 min at 60°C. The reaction leads to the mono-derivatives. An aliquot of the derivatized solution is used for GC/MSD measurement.


Caution: Only diluted samples are injected in order to avoid column stationary phase deterioration.

Results


Underivatized Barbiturates can be measured without problems. The derivatization improves sensitivity and also increases molecular mass (+181amu) which is advantageous especially in SIM mode. A drastic improvement is noticed in ECNI mode. The signal/noise ratio for 1pg/µl analyte concentration in ECNI mode exceeds 200:1. Sensitivity for the derivatives in PCI with NH₃ is very low.




El-Spectrum, Amobarbital, underivatized, m/z 226; M⁺




El-Spectrum, Amobarbital, PFBB-derivative, m/z 406; M⁺




PCI/CH₄-Spectrum, Amobarbital, underivatized: *m/z* 227, 255, 267; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




PCI/CH₄-Spectrum, Amobarbital, PFBB-derivative: *m/z* 407, 435, 447; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




ECNI/CH₄-Spectrum, Amobarbital, PFBB-derivative: *m/z* 405; [M-H]⁻




EI-Spectrum, Barbitol, underderivatised, m/z 184; M^+




EI-Spectrum, Barbitol, PFBB-derivative, m/z 364; M^+




PCI/ CH_4 -Spectrum, Barbitol, underderivatised: m/z 185, 213, 225; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




PCI/CH₄-Spectrum, Barbitol, PFBB-derivative: *m/z* 365, 393, 405; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




ECNI/CH₄-Spectrum, Barbitol, PFBB-derivative: *m/z* 363; [M-H]⁻




EI-Spectrum, Butethal, underivatized, *m/z* 212; M⁺




EI-Spectrum, Butethal, PFBB-derivative, m/z 392; M^+




PCI/ CH_4 -Spectrum, Butethal, underderivatised: m/z 213, 241, 253; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




PCI/ CH_4 -Spectrum, Butethal, PFBB-derivative: m/z 393, 421, 433; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




ECNI/CH₄-Spectrum, Butethal, PFBB-derivative: m/z 391; [M-H]⁻




El-Spectrum, Pentobarbital, underderivatised, m/z 226; M⁺




El-Spectrum, Pentobarbital, PFBB-derivative, m/z 406; M⁺




PCI/CH₄-Spectrum, Pentobarbital, underderivatised: *m/z* 227, 255, 267; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




PCI/CH₄-Spectrum, Pentobarbital, PFBB-derivative, *m/z* 407, 435, 447; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




ECNI/CH₄-Spectrum, Pentobarbital, PFBB-derivative: *m/z* 405; [M-H]⁻




EI-Spectrum, Secobarbital, underivatised, m/z 238; M⁺




EI-Spectrum, Secobarbital, PFBB-derivative, m/z 418; M⁺



PCI/CH₄-Spectrum, Secobarbital, underivatised: m/z 239, 267, 279; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/CH₄-Spectrum, Secobarbital, PFBB-derivative: m/z 419, 447, 459; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




ECNI/CH₄-Spectrum, Secobarbital, PFBB-derivative: m/z 417; $[M-H]^-$

ECNI SIM, 5 Barbiturates

Analytes in elution order	RT (min)	Nominal ions (m/z)	Relative Signal/Noise Ratio
Barbital	14.52	363	1
Butethal	15.37	391	1.03
Amobarbital	15.63	405	1.00
Pentobarbital	15.93	405	1.15
Secobarbital	16.13	417	0.72

Table 1



ECNI SIM, 5 Barbiturates, 1 μ g/ μ l each, see Table 1

Benzodiazepines

Alprazolam CAS-Nr. 28981-97-7

Molecular formula: C₁₇H₁₃ClN₄

Bromazepam CAS-Nr. 1812-30-2

Molecular formula: C₁₄H₁₀BrN₃O

Diazepam CAS-Nr. 439-14-5

Molecular formula: C₁₆H₁₃ClN₂O

Flunitrazepam CAS-Nr. 1622-62-4

Molecular formula: C₁₆H₁₂FN₃O₃

Triazolam CAS-Nr. 28911-01-5

Molecular formula: C₁₇H₁₂C₁₂N₄

GC-Parameters

Column: HP-5ms

Agilent Part Nr.19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection:

On Column: 100°C/40°C

Oven Temp.Program

Scan: 100°C (0.3min)

SIM: 40°C (0.3min) -

25°C/min to 300°C (6min)

Flunitrazepam

Scan/SIM: 100°C (0.3min)

MS-Parameter

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Remarks

Injection

“Active” surfaces in the inlet system may causes discrimination. For these measurements, on-column injection technique with a fused silica needle syringe was applied to allow a better comparison.

ECNI Parameters

The commonly applied parameters for improving ECNI measurements (like Emission Current, Flow, NH₃ Buffer Gas) showed no positive effects.


Results

Sensitivity in ECNI Scan mode was measured using concentrations in the range of 0.3ng/µl to 1ng/µl and is in that range. SIM mode measurements were made at about 1000 times lower; 3pg/µl to 10pg/µl. Flunitrazepam showed the highest relative response (at < 200fg/µl).


References

“Application of Electron Capture Negative Chemical Ionization for the detection of a Date Rape Drug”


A. Negrusz, Ch. Moore, H. Prest
Agilent Pub. Nr. 5968-4364E




EI-Spectrum, Alprazolam: m/z 308; M⁺




ECNI-Spectrum, Alprazolam: m/z 308; M⁻




El-Spectrum, Bromazepam: m/z 315; M^+




ECNI-Spectrum, Bromazepam: m/z 315; M^+




El-Spectrum, Diazepam: m/z 284; M^+




ECNI-Spectrum, Diazepam: m/z 284; M^+




EI-Spectrum, Flunitrazepam: m/z 312; M^+



ECNI-Spectrum, Flunitrazepam: m/z 313; M^+




El-Spectrum, Triazolam: m/z 342; M^+



ECNI-Spectrum, Triazolam: m/z 342; M^-

ECNI SIM




ECNI SIM, Alprazolam


ECNI SIM

Analytes	Retention Time (min)	Ions (m/z)	Signal/Noise Relative to Diazepam (TIC)
Alprazolam	11.49	308, 310	0.6
Bromazepam	9.60	315, 317	0.08
Diazepam	8.77	284, 286	1
Flunitrazepam	9.49	313	16
Triazolam	12.32	306, 308	0.7


Table: Sensitivity Benzodiazepines




Bromazepam



Diazepam



Flunitrazepam



Triazolam

Benzoylecgonine

CAS-Nr. 519-09-5

Molecular formula: C₁₆H₁₉NO₄

GC-Parameters

Column: HP-5ms

Agilent Part Nr.19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to

300°C (5min)

MS-Parameter

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Remarks

Derivatization

a) Trimethyl silylation –TMS– with MSTFA (Reagent: Fluka 69479)

b) Reaction with Pentafluoropropionic Acid Anhydride (PFPA) (Reagent: Fluka 77292)

a) The standard solution (SIGMA B 8900), concentration 100ng/µl, diluted in ethyl acetate, was evaporated with a gentle nitrogen flow. To the residue, 50µl reagent is added and the reaction mixture is incubated for 20min at 60°C.

Gentle evaporation with nitrogen is repeated and the residue redissolved in ethyl acetate.


b) Procedure as described above. After the first evaporation the residue, 80µl of the PFPA reagent and 20µl of hexafluoro-isopropanol (FLUKA 52517) are added and the reaction mixture is incubated for 30 min at 70°C. Then evaporation, dilution and GC/MSD analysis.

Results


Derivatization is recommended.

In PCI/NH₃ Mode, TMS Derivative, the degree of fragmentation is related to sample concentration for this analyte.


Sensitivity is directly related to the derivatization and to the applied reagent gas, see table. In SIM mode, analyte concentration of 1pg/µl is measured with signal/noise ratio of approximately 10/1.




EI-Spectrum, Benzoylecgonine, underivatized: m/z 289; M⁺




EI-Spectrum, Benzoylecgonine, TMS Derivative: m/z 361; M⁺




El-Spectrum, Benzoyllecgonine, PFPA Derivative, -O-CH-(CF₃)₂: m/z 439; M⁺




PCI/CH₄-Spectrum, Benzoyllecgonine, 50ng/μl, TMS Derivative: m/z 362, 390, 402; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




PCI/CH₄-Spectrum, Benzoyllecgonine, 10ng/μl, PFPA Derivative, -O-CH-(CF₃)₂: m/z 440, 468, 480; [M+H]⁺, [M+C₂H₅]⁺, [M + C₃H₅]⁺



PCI/NH₃-Spectrum, Benzoylecgonine, 10ng/μl, TMS Derivative: m/z 362; [M + H]⁺




PCI/NH₃-Spectrum, Benzoylecgonine, 50ng/μl, TMS Derivative: m/z 362; [M+H]⁺



PCI/NH₃-Spectrum, Benzoylecgonine, 10ng/μl, PFP Derivative, -O-CH-(CF₃)₂: m/z 440; [M+H]⁺

PCI/NH₃ – SIM Mode



**Benzoylcgonine, PFPA Derivative,
Retention Time: 8.84min
1pg/μl, SIM Ions: 318, 440 m/z**

Derivative	Ion - Mode	Peak to Peak Signal:Noise Ratio
TMS	PCI/CH ₄	> 10 : 1
TMS	PCI/NH ₃	> 35 : 1
PFPA	PCI/CH ₄	> 100 : 1
PFPA	PCI/NH ₃	> 200 : 1

**Table : Benzoylcgonine, Sensitivity
(S/N), Scan Acquisition, 10ng/μl each**

Chloramphenicol

CAS-Nr. 56-75-7

Molecular Formula: $C_{11}H_{12}Cl_2N_2O_5$

CAS-Nr. O,O-TMS Derivative:

21196-84-9

GC-Parameter

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25 μ m

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Pulsed splitless

Oven Temp. Program

70°C (1min) - 30°C/min to 150°C

15°C/min to 300°C

MS-Parameter

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1.75ml/min (35)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

SIM Tune+ 600V

Remarks

Derivatization

Trimethyl silylation with

HMDS/TMCS at 2/1 in Pyridine

(Reagent: Fluka 85431)

The standard solution (SIGMA

C 0378), concentration 20ng/ μ l,

diluted in ethyl acetate, is gently

evaporated under dry nitrogen.

To the residue, 50 μ l reagent is

added and the reaction mixture

is incubated for 2 min at 50°C.

Evaporation to dryness is repeated

and the residue diluted in hexane.

The solution is ready for injection.

Results

Analyzing this compound as

O,O TMS derivative is preferred

over the parent.

PCI Scan sensitivity is on the

order of 20ng/ μ l with methane

however ammonia produces

5 times greater sensitivity.

The greatest sensitivity is


achieved in ECNI SIM mode

which allows detection at 0.1 pg/ μ l.


Calibration curves show good

linearity over the concentration


range of 0.1pg/ μ l to 100pg/ μ l.




EI Spectrum, Chloramphenicol, O,O-TMS Derivative: m/z 466; M⁺



PCI/CH₄ Spectrum, Chloramphenicol, O,O-TMS Derivative: m/z 467, 495, 507; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




PCI/NH₃ Spectrum, Chloroamphenicol, O,O-TMS Derivative: m/z 467, 484; [M+H]⁺, [M+NH₄]⁺




ECNI/CH₄ Spectrum, Chloroamphenicol, O,O-TMS Derivative: m/z 466; M⁻



SIM Mode (Ions: 465.9, 467.9, 469.9 m/z)



ECNI/CH₄ Chloroamphenicol, O,O-TMS Derivative, 1pg/μl



ECNI/CH₄ Chloroamphenicol, O,O-TMS Derivative, 0.1pg/μl



Calibration Curves: ECNI/CH₄ Chloroamphenicol, 0,0-TMS Derivative
 left: 0.1pg/μl to 100pg/μl, right: 0.1pg/μl to 1pg/μl

Chlorphenoxamine

CAS-Nr. 77-38-3

Molecular Formula: C₁₈H₂₂ClNO

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI - SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 15.04min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 500/1


PCI/NH₃ Scan: >100/1



EI-Spectrum, Chlorphenoxamine: m/z 303; M^+



PCI/CH₄-Spectrum, Chlorphenoxamine: m/z 304, 332, 344; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$



PCI/NH₃-Spectrum, Chlorphenoxamine: m/z 304; $[M+H]^+$

Chlorprothixene

CAS-Nr. 113-59-7

Molecular Formula: C₁₈H₁₈ClNS

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 10.56min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 25/1


PCI/NH₃ Scan: > 60/1



EI-Spectrum, Chlorprothixene: m/z 315; M⁺



PCI/CH₄-Spectrum, Chlorprothixene: m/z 316, 344, 356; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃-Spectrum, Chlorprothixene: m/z 316 : [M+H]⁺

Cimaterol

CAS-Nr. 54239-37-1

Molecular Formula: C₁₂H₁₇N₃O

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to

280°C (5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltages: Tune + 400V

Remarks

Derivatizations: TMS and PFPA

a) Trimethyl silylation (TMS)

with BSTFA/TMCS

(Reagent: Fluka 15238)

The standard solution (Boehringer Ingelheim), concentration 1mg/ml, diluted in methanol, is evaporated to dryness with a gentle nitrogen flow. To the residue, 50µl reagent and 125 µl pyridine is added and the reaction mixture is incubated for 30min at 60°C. Evaporation is repeated and the residue is redissolved in chloroform. The solution is ready for injection.

b) Reaction with Pentafluoropropionic Acid Anhydride (PFPA)

(Reagent: Fluka 77292)

The standard is treated as in a)

above. After the first evaporation

the residue is treated with 80µl

of the PFPA reagent, 20µl of

Hexafluoroisopropanol (Fluka

52517) is added and the reaction

mixture is incubated for 30min at

70°C, followed by evaporation,

dilution and GC/MSD analysis.

Results

The EI spectra of underivatized and derivatized analytes show low intensity for the molecular ion.

Both the TMS and PFPA

derivatization reactions form

the di-derivatives. The TMS

derivative response is higher

comparing to the PFPA derivative.

PCI/NH₃ response shows

increased sensitivity by factor 2.6.

Both modes present the

characteristic PCI adduct ion


information. Sensitivity for the

TMS derivative in PCI/NH₃ SIM


mode results in a signal-to-noise

ratio of >20:1 at an analyte


concentration of 50pg/µl.




EI-Spectrum, Cimaterol, underivatized: m/z 219; M⁺




El-Spectrum, Cimaterol, TMS Derivative: m/z 363; M^+




El-Spectrum, Cimaterol, PFPA Derivative: m/z 511; M^+




PCI/ CH_4 -Spectrum, Cimaterol, underderivatised: m/z 220, 248, 260; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$; m/z 202, 230, 242; $[M-OH]^+$, $[M+C_2H_5-OH_2]^+$, $[M+C_3H_5-OH_2]^+$




PCI/CH₄-Spectrum, Cimaterol, TMS Derivative: *m/z* 364, 392, 404; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




PCI/NH₃-Spectrum, Cimaterol, TMS Derivative: *m/z* 364 [M+H]⁺



PCI/CH₄-Spectrum, Cimaterol, PFP Derivative: *m/z* 494, 522, 534; [M-OH]⁺, [M+C₂H₅-OH₂]⁺, [M+C₃H₅-OH₂]⁺




PCI/NH₃-Spectrum, Cimaterol, PFPA Derivative: m/z 529; [M+NH₄]⁺



ECNI/CH₄-Spectrum, Cimaterol, PFPA Derivative: m/z 491; [M-H]⁻

PCI/NH₃ – SIM Mode



Cimaterol, TMS Derivative, 50pg, Retention Time: 8.98min
Ions: 274, 364 m/z ; Signal/Noise: > 20/1

Clenbuterol

CAS-Nr. 37148-27-9

Molecular Formula: C₁₂H₁₈C₁₂N₂O

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Split & Pulsed splitless, 250°C

Oven Temp. Program

Split: 100°C (0.3min) –

25°C/min to 280°C

Pulsed splitless: 80°C (1min) –

25°C/min to 265°C (1min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN/SIM

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN/SIM

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Remarks

Derivatization

Trimethyl silylation with BSTFA/TMCS (Reagent: Fluka 15238)

Derivatization conditions (e.g, incubation temperature and time) have some impact on the product formation. Applying 60°C for 30min generates both mono- and


di-derivatives at a ratio of 1:0.75. The standard solution (SIGMA C 5423), concentration 1mg/ml, in methanol, is evaporated with a gentle nitrogen flow. To the residue a mixture of Pyridine/TMS-Reagent (2.5/1) is added and incubation is done for 30min. The excess reagent is evaporated and the residue redissolved in chloroform. The solution is ready for injection.

Results


The degree of fragmentation is related to the sample concentration for this analyte. The response in PCI/NH₃ mode is higher comparing to PCI/CH₄ mode. Weak response is observed in ECNI/CH₄ mode.

References


“Analysis of Clenbuterol by GC-MS...”, F. David, RIC, Agilent Pub. Nr. 5962-9427E
“Clenbuterol and Norandrosterone...”, B. Wüst, Agilent Pub. Nr. 5980-0908E




EI-Spectrum, Clenbuterol, underivatized: m/z 276; M⁺




EI-Spectrum, Clenbuterol, BSTFA-mono-derivative: m/z 348; M⁺




El-Spectrum, Clenbuterol, BSTFA-di-derivative: m/z 420; M^+




PCI/ CH_4 -Spectrum, Clenbuterol, underivatized: m/z 277, 305, 317; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




PCI/ CH_4 -Spectrum, Clenbuterol, BSTFA-mono-derivative, approx. 30ng: m/z 349, 377, 389; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




PCI/CH₄-Spectrum, Clenbuterol, BSTFA-mono-derivative, approx. 1ng; m/z 349, 377, 389; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




PCI/CH₄-Spectrum, Clenbuterol, BSTFA-di-derivative, approx. 40ng; m/z 421, 449, 461; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




PCI/CH₄-Spectrum, Clenbuterol, BSTFA-di-derivative, approx. 1ng; m/z 421, 449, 461; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$




PCI/NH₃-Spectrum, Clenbuterol, underivatised: m/z 277; [M+H]⁺




PCI/NH₃-Spectrum, Clenbuterol, mono-derivative, approx. 30ng: m/z 349; [M+H]⁺




PCI/NH₃-Spectrum, Clenbuterol, mono-derivative, approx. 1ng: m/z 349; [M+H]⁺




PCI/NH₃-Spectrum, Clenbuterol, approx. 40ng, di-derivative: m/z 421; [M+H]⁺



PCI/NH₃-Spectrum, Clenbuterol, di-derivative, approx. 1ng: m/z 421; [M+H]⁺




ECNI/CH₄-Spectrum, Clenbuterol, approx. 1ng, underivatized: m/z 276; M⁻




ECNI/CH₄-Spectrum, Clenbuterol, approx. 1ng, di-derivative: m/z 420; M⁻


SIM Mode



PCI/CH₄, Clenbuterol, approx. 15pg, mono & di-derivative: S/N ≈ 70/1 & ≈ 5/1



PCI/NH₃, Clenbuterol, approx. 1.5pg, mono & di-derivative: S/N ≈ 55/1 & ≈ 8/1



ECNI/CH₄, Clenbuterol, approx. 15 pg, di-derivative: S/N ≈ 30/1

Cocaine

CAS-Nr. 50-36-2

Molecular Formula: C₁₇H₂₁NO₄

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp.Program

70°C (1min) – 25°C/min to

300°C (5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia


Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File


EM Voltage: Tune + 400V

Results


In PCI mode, SCAN acquisition, ammonia is the preferred reagent gas due to higher response compared to methane. In SIM mode there were no significant differences in responses observed. In PCI/CH₄ SIM mode signal/noise ration of >20/1 was measured for 10pg/µl. ECNI showed no relevant spectrum.



EI-Spectrum, Cocaine: m/z 303: M⁺




PCI/CH₄-Spectrum, Cocaine: m/z 304, 332, 344; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃-Spectrum, Cocaine: m/z 304; [M+H]⁺

PCI/CH₄ – SIM Mode



Cocaine, Retention Time: 9.50min, 10pg/μl: Ions 182, 304 m/z ;
S/N > 25/1

Codeine

CAS-Nr. 76-57-3

Molecular Formula: C₁₈H₂₁NO₃

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to 300°C

(5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Remarks

Derivatization

Reaction with Pentafluoropropionic Acid Anhydride (PFPA)


(Reagent: Fluka 77292)

The standard solution (SIGMA


C 1653), concentration 100ng/µl in ethyl acetate, is evaporated with a gentle flow of nitrogen. To the residue 80µl of the PFPA reagent and 20µl of Hexafluoroisopropanol (Fluka 52517) is added and the reaction mixture is incubated for 30min at 70°C. Evaporation is repeated and the residue redissolved in ethyl acetate. The solution is ready for injection and analysis.

Results


Even the underivatized analyte is measured without chromatographic discrimination. In PCI mode, the response of the derivatized analyte was 2-times higher using ammonia reagent gas than methane. ECNI/CH₄ SIM measurements were highly sensitive for the derivatized analyte; 5pg/µl resulted in signal/noise ratio of >125/1. The acetylated analyte showed no significant ECNI spectrum.




EI-Spectrum, Codeine, underivatized: m/z 299 ; M⁺




EI-Spectrum, Codeine, PFPA Derivative: m/z 445 ; M⁺



PCI/CH₄-Spectrum, Codeine, PFPA Derivative: m/z 446, 474, 486; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺




PCI/NH₃-Spectrum, Codeine, PFPA Derivative: m/z 446, 463; [M+H]⁺, [M+NH₄]⁺



ECNI/CH₄-Spectrum, Codeine, PFPA Derivative: m/z 425; [M-H]⁻

ECNI/CH₄ SIM



**Codeine, PFPA Derivative, 5pg, Retention Time: 10.08min
Ion: 425 m/z, Signal/Noise \approx 150/1**

Dimethindene

CAS-Nr. 5636-83-9

Molecular Formula: C₂₀H₂₄N₂

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 9.59min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 5/1


PCI/NH₃ Scan: > 10/1



EI-Spectrum, Dimethindene: *m/z* 292; M⁺



PCI/CH₄-Spectrum, Dimethindene: *m/z* 293, 321, 333; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃-Spectrum, Dimethindene: *m/z* 293; [M+H]⁺

Diphenhydramine

CAS-Nr. 58-73-1

Molecular Formula: C₁₇H₂₁NO

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 12.81min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: >500/1


PCI/NH₃ Scan: >80/1



EI-Spectrum, Diphenhydramine: m/z 255; M⁺



PCI/CH₄-Spectrum, Diphenhydramine: m/z 256, 284, 296; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃-Spectrum, Diphenhydramine: m/z 256; [M+H]⁺

Lidocaine

CAS-Nr. 137-58-6

Molecular Formula: C₁₄H₂₂N₂O

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperature:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 12.95min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: >350/1


PCI/NH₃ Scan: >80/1



EI-Spectrum, Lidocaine: *m/z* 234; M⁺



PCI/CH₄-Spectrum, Lidocaine: *m/z* 235, 263, 275; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃-Spectrum, Lidocaine: *m/z* 235; [M+H]⁺

Mabuterol

CAS-Nr. 56341-08-3

Molecular Formula: C₁₃H₁₈ClF₃N₂O

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min

to 280°C (5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN/SIM

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Remarks

Derivatization

a) Trimethyl silylation with

BSTFA/TMCS

(Reagent: Fluka 15238)

b) Reaction with Pentafluoropro-

panionic Acid Anhydride – PFPA –

(Reagent: Fluka 77292)


a) 100µl of the hydrochloride standard (Boehringer Ingelheim), concentration 1.6mg/ml, is dissolved in methanol and evaporated with a gentle flow of nitrogen. To the residue, 50µl of derivatization reagent and 125µl of dry pyridine are added and the reaction mixture is incubated for 30min at 60°C. Gentle evaporation with nitrogen is repeated and the residue dissolved in chloroform.

b) Procedure as described above. After the first evaporation to the residue, 80µl of the derivatization reagent and 20µl of hexafluoro isopropanol (Fluka 52517) are added and the reaction mixture is incubated for 30min at 70°C. Then evaporation, dilution and GC/MSD analysis.


Results

The EI spectra of the underivatized and derivatized analytes show low intensity for the molecular ion. The TMS derivatization react to form a mono-derivative and the PFPA derivatization forms a di-derivative. The PFPA spectra show M-18 fragmentation. In PCI/CH₄ mode the derivatives generate the molecular ions and show the characteristic adduct ions. The measurements of the TMS derivatives in PCI/NH₃ Scan mode are by factor 6 more sensitive than the PCI/CH₄ measurements. The PFPA derivative shows less response than the TMS derivative in PCI mode.


An improvement is noticed in ECNI mode for the PFPA derivative. In SIM mode, the signal/noise ratio is approximately 40:1 for an analyte concentration of 200fg/µl.




El-Spectrum, Mabuterol, underivatized, m/z 310; M⁺




El-Spectrum, Mabuterol, TMS mono-derivative, m/z 382; M⁺




El-Spectrum, Mabuterol, PFPA di-derivative, m/z 602; M⁺




PCI/CH₄-Spectrum, Mabuterol, underivatised, m/z 311, 339, 351; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺




PCI/CH₄-Spectrum, Mabuterol, TMS mono-derivative, *m/z* 383, 411, 423; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺




PCI/NH₃-Spectrum, Mabuterol, TMS mono-derivative, *m/z* 383; [M + H]⁺



PCI/CH₄-Spectrum, Mabuterol, PFPA di-derivative, *m/z* 585, 613, 625; [M-18 + H]⁺, [M-18 + C₂H₅]⁺, [M-18 + C₃H₅]⁺




PCI/NH₃-Spectrum, Mabuterol, PFPA di-derivative, m/z 602; M⁺



ECNI/CH₄-Spectrum, Mabuterol, PFPA di-derivative, molecular mass = 602 u

ECNI/CH₄ – SIM Mode




Mabuterol, PFPA derivative, 200fg,
Retention Time: 8.24min
Ions: m/z 507/509, Signal/Noise: 40/1

MDA

3,4-Methylenedioxyamphetamine

CAS-Nr. 4764-17-4

Molecular Formula: C₁₀H₁₃NO₂



GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to 300°C (5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode:

ECNI/CH₄/NH₃ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Buffer Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Remarks

Derivatization

- Trifluoroacetylation (TFA) with MBTFA (Reagent: Fluka 65943)
- Reaction with Pentafluoropropionic Acid Anhydride – PFPA – (Reagent: Fluka 77292)

a) 100µl of the standard (SIGMA M 3272), 100ng/µl MDA concentration, is dissolved in ethyl acetate and evaporated with a gentle flow of nitrogen. To the residue 50µl derivatization


reagent is added and the solution is incubated for 30min at 80°C.

Gentle evaporation with nitrogen is repeated and the residue dissolved in ethyl acetate.


b) Procedure as described above. After the first evaporation to the residue, 80µl of the derivatization reagent and 20µl of hexafluoroisopropanol (Fluka 52517) are added and the reaction mixture is incubated for 30min at 70°C. Then evaporation, dilution and GC/MSD analysis.

Results


Derivatization is recommended for this analyte. In PCI/NH₃ mode the analyte spectrum shows the NH₄ adduct ion as base peak. Signal is a factor 1.5 more sensitive compared to PCI with CH₄. In ECNI mode the TFA derivative response is relatively low. PFPA derivatization is the reaction of choice for both PCI and ECNI modes. Fragmentation and sensitivity are related to the choice of derivatizing reagent and buffer gas. In ECNI/NH₃ SIM mode the signal:noise ratio is approximately 65:1 for the 5pg PFPA derivative.




El-Spectrum, MDA, TFA derivative, m/z 275; M⁺




PCI/CH₄-Spectrum, MDA, TFA derivative, *m/z* 276, 304, 316; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺




PCI/NH₃-Spectrum, MDA, TFA derivative, *m/z* 293; [M + NH₄]⁺




EI-Spectrum, MDA, PFPA derivative, *m/z* 325; M⁺




PCI/CH₄-Spectrum, MDA, PFPA derivative, *m/z* 326, 354, 366; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺



PCI/NH₃-Spectrum, MDA, PFPA derivative, *m/z* 343; [M + NH₄]⁺




ECNI/CH₄-Spectrum, MDA, PFPA derivative, *m/z* 305, 324; [M - HF]⁻, [M - H]⁻



ECNI/NH₃-Spectrum, MDA, PFPA derivative, m/z 305, 324; [M-HF]⁻, [M-H]⁻

ECNI/NH₃ – SIM Mode



MDA, PFPA derivative, 5pg,
Retention Time: 7.12min
Ions: m/z 305, 324; Signal: Noise ≈ 65/1

Mepivacaine

(Carbocaine)

CAS-Nr. 96-88-8

Molecular Formula: $C_{15}H_{22}N_2O$

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25 μ m

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperature:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 14.89min

Analyte Concentration: 4ng/ μ l


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 680/1


PCI/NH₃ Scan: > 70/1



EI-Spectrum, Mepivacaine: m/z 246; M⁺



PCI/CH₄-Spectrum, Mepivacaine: m/z 247, 275, 287; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃-Spectrum, Mepivacaine: m/z 247; [M+H]⁺

Methadone

CAS-Nr. 76-99-3

Molecular Formula: C₂₁H₂₇NO

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 15.54min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 300/1


PCI/NH₃ Scan: > 30/1



EI-Spectrum, Methadone: m/z 309; M⁺



PCI/CH₄-Spectrum, Methadone: m/z 310, 338, 350; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃-Spectrum, Methadone: m/z 310; [M+H]⁺

Morphine

CAS-Nr. 57-27-2

Molecular Formula: C₁₇H₁₉NO₃

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to 300°C

(5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Remarks

Derivatization

Reaction with Pentafluoropropionic Acid Anhydride (PFPA)

(Reagent: Fluka 77292)

100µl of the standard (SIGMA M 9524), concentration 100ng/µl, dissolved in ethyl acetate is


evaporated with a gentle flow of nitrogen. To the residue 80µl derivatization reagent and 20µl hexafluoroisopropanol are added and the mixture is incubated for 30min at 70°C. Gentle evaporation with nitrogen is repeated and the residue redissolved in ethyl acetate. The solution is ready GC/MSD analysis.

Results


Derivatization is recommended.

The trifluoroacetylation (TFA) with MBTFA leads in EI mode to a moderately intense molecular ion. In ECNI mode the TFA analyte spectrum shows no distinctive results for the derivative.


The PFPA derivatization is suitable for both PCI/NH₃ and for ECNI/CH₄ measurements. The signal-to-noise ratio for 1pg analyte PFPA derivative in ECNI/CH₄ mode is approximately 30:1.




El-Spectrum, Morphine, underivatized, m/z 285; M⁺




El-Spectrum, Morphine, Trifluoroacetyl derivative, m/z 477; M⁺




El-Spectrum, Morphine, PFPA derivative, m/z 577; M⁺



PCI/CH₄-Spectrum, Morphine, PFPA derivative, m/z 578, 606, 618; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺




PCI/NH₃-Spectrum, Morphine, PFPA derivative, m/z 578; [M + H]⁺



ECNI/CH₄-Spectrum, Morphine, PFPA derivative, *m/z* 537, 557, 577; [M – 2(HF)]⁻; [M-HF]⁻; [M]⁻

ECNI/CH₄- SIM Mode



Morphine, PFPA derivative, Retention Time: 9.89min
1pg/μl, Ions: *m/z* 537, 557; S/N ≈ 30/1

Acquisition Mode	Analyte conc.	Approximate Signal/Noise Ratio
EI-Scan	10ng/μl	250/1
PCI/CH ₄ -Scan	10ng/μl	70/1
PCI/NH ₃ -Scan	10ng/μl	150/1
ECNI/CH ₄ -Scan	10ng/μl	110/1
ECNI/CH ₄ -SIM	1pg/μl	30/1

Table: Morphine, PFPA derivative, Sensitivity (S/N), EI/PCI/ECNI, Scan/SIM

Nalorphine

CAS-Nr. 62-67-9

Molecular Formula: C₁₉H₂₁NO₃

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to 300°C

(5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Remarks

Derivatization

Reaction with Pentafluoropropionic

Acid Anhydride (PFPA)

(Reagent: Fluka 77292)

100µl of the standard (SIGMA

N 0762), concentration 100ng/µl,

dissolved in ethyl acetate is


evaporated with a gentle flow of nitrogen. To the residue, 80µl of derivatization reagent and 20µl hexafluoroisopropanol are added and the mixture is incubated for 30min at 70°C. Gentle evaporation with nitrogen is repeated and the residue redissolved in ethyl acetate. The solution is ready GC/MSD analysis.

Results


Derivatization is recommended.

The trifluoroacetylation (TFA) with MBTFA leads in EI mode to a moderately intense molecular ion. In ECNI mode the TFA analyte spectrum exhibits no distinctive features.


The PFPA derivatization is suitable for PCI/NH₃ and for ECNI/CH₄ measurements. The signal/noise ratio for 1pg analyte PFPA derivative in ECNI/CH₄ mode is approximately 15:1.




El-Spectrum, Nalorphine, underivatized, m/z 311; M⁺




El-Spectrum, Nalorphine, Trifluoroacetyl derivative, m/z 503; M⁺




EI-Spectrum, Nalorphine, PFPA derivative, m/z 603; M⁺



PCI/CH₄-Spectrum, Nalorphine, PFPA derivative, m/z 604, 632, 644; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺




PCI/NH₃-Spectrum, Nalorphine, PFPA derivative, m/z 604; [M + H]⁺



ECNI/CH₄-Spectrum, Nalorphine, PFPA derivative, *m/z* 563, 583, 603; [M – 2(HF)]⁻, [M – HF]⁻, M⁻

ECNI/CH₄- SIM Mode



Nalorphine, PFPA derivative, Retention Time: 10.30min
1pg/μL, Ions: *m/z* 563, 583; S/N: 16/1

Acquisition Mode	Analyte conc.	Approximate Signal/Noise Ratio
EI-Scan	10ng/μl	240/1
PCI/CH ₄ -Scan	10ng/μl	75/1
PCI/NH ₃ -Scan	10ng/μl	110/1
ECNI/CH ₄ -Scan	10ng/μl	200/1
ECNI/CH ₄ -SIM	1pg/μl	15/1

Table: Nalorphine, PFPA derivative, Sensitivity (S/N), EI/PCI/ECNI, Scan/SIM

Nitroimidazoles

Dimetridazole

CAS-Nr. 551-92-8

Molecular Formula: C₅H₇N₃O₂

Ronidazole

CAS-Nr. 7681-76-7

Molecular Formula: C₆H₈N₄O₄

Metronidazole

CAS-Nr. 443-48-1

Molecular Formula: C₆H₉N₃O₃

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

60°C (1min) – 25°C/min to 270°C

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Buffer Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage Scan: Tune + 300V

EM Voltage SIM: Tune + 400V

Remarks

Derivatization

Ronidazole (Reagent: SIGMA

R 7635) and Metronidazole

(Reagent: SIGMA M 1547)

can be silylated with MSTFA

(Reagent: Fluka 69479)

Each standard solution,

concentration 1ng/µl in ethyl

acetate, is evaporated with a gentle

nitrogen flow. To the residue, 50µl

of the reagent is added and

incubated at 15min at 60°C.

The derivatized solutions are used for GC/MSD measurements.

Results


The spectrum of Ronidazole shows no molecular ions either in EI or in ECNI mode. Spectral base peaks of the analytes are related to the fragmentation at the carbamide ester positions.

The degree of fragmentation is affected by choice of the buffer gas. Both methane and ammonia show different responses.


Considering absolute signal intensity, ammonia is the preferred buffer gas. However the signal/noise ratios indicate that both gases perform comparably. The relative S/N ratios in SIM mode for these nitroimidazoles are approximately:

Dimetridazole : Metronidazole :


Ronidazole = 100 : 10 : 1.




Dimetridazole




Ronidazole




Metronidazole




EI-Spectrum, Dimetridazole, m/z 141; M⁺




ECNI/CH₄-Spectrum, Dimetridazole, m/z 141; M⁺




El-Spectrum, Ronidazole, TMS derivative, molecular mass = 272 u




ECNI/CH₄-Spectrum, Ronidazole, TMS derivative, molecular mass = 272 u



ECNI/NH₃-Spectrum, Ronidazole, TMS derivative, molecular mass = 272 u




EI-Spectrum, Metronidazole, TMS derivative, m/z 243; M⁺




ECNI/CH₄-Spectrum, Metronidazole, TMS derivative, m/z 243; M⁺

Scan Mode, ECNI, 1ng each




TIC of Dimetridazole(1), Ronidazole(2), Metronidazole(3)
Buffer Gas: Methane, bottom; Ammonia, top


SIM Mode, ECNI/CH₄, 1 pg each



Ronidazole, S/N \approx 6/1
Ions m/z 139/229



Metronidazole, S/N \approx 60/1
Ion m/z 243



Dimetridazole, S/N \approx 600/1
Ion m/z 141

Orphenadrine

CAS-Nr. 83-98-7

Molecular Formula: C₁₈H₂₃NO

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperature:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 13.54min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 300/1


PCI/NH₃ Scan: > 50/1



EI-Spectrum, Orphenadrine: m/z 269; M⁺



PCI/CH₄-Spectrum, Orphenadrine: m/z 270, 298, 310; [M+H]⁺, [M+C₂H₅]⁺, [M+C₃H₅]⁺



PCI/NH₃-Spectrum, Orphenadrine: m/z 270; [M+H]⁺

Phenylbutazone

CAS-Nr. 50-33-9

Molecular Formula: C₁₉H₂₀N₂O₂

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to 300°C

(2min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:


Source 150°C, Quad 106°C

EM Voltage: Tune + 400V


Remarks

Results

In EI mode, the analyte shows a moderately intense molecular ion. ECNI/CH₄ mode generates the molecular anion and SIM measurements in this mode result in an approximate signal/noise ratio of 15:1 for 1pg of analyte.




EI-Spectrum, Phenylbutazone, m/z 308; M⁺



ECNI/CH₄-Spectrum, Phenylbutazone, m/z 308; M⁻

ECNI/CH₄ – SIM Mode



Phenylbutazone, 1pg, Retention Time: 10.32min
Ion: m/z 308, Signal/Noise ≈ 7/1

Promethazine

CAS-Nr. 60-87-7

Molecular Formula: $C_{17}H_{20}N_2S$

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25 μ m

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 16.79min

Analyte Concentration: 4ng/ μ l


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 200/1


PCI/NH₃ Scan: > 90/1



EI-Spectrum, Promethazine: m/z 284 ; M^+



PCI/CH₄-Spectrum, Promethazine: m/z 285, 313, 325; $[M+H]^+$, $[M+C_2H_5]^+$, $[M+C_3H_5]^+$



PCI/NH₃-Spectrum, Promethazine: m/z 285; $[M+H]^+$

Propionylpromazine

Combelen

CAS-Nr. 3568-24-9

Molecular Formula: C₂₀H₂₄ON₂S

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Split, 250°C

Oven Temp. Program

120°C (0.3min) – 20°C/min to

300°C (4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Results


Analyte Retention Time: 9.25min

Analyte Concentration: 4ng/µl


Signal/Noise Ratios for TIC

PCI/CH₄ Scan: > 8/1


PCI/NH₃ Scan: > 16/1



EI-Spectrum, Propionylpromazine: *m/z* 340; M⁺



PCI/CH₄-Spectrum, Propionylpromazine: *m/z* 327, 355, 367; [M+H-CH₂]⁺, [M+C₂H₅-CH₂]⁺, [M+C₃H₅-CH₂]⁺



PCI/NH₃-Spectrum, Propionylpromazine: *m/z* 341; [M+H]⁺

Ractopamine

CAS-Nr. 99095-19-9 (HCl)

Molecular Formula: C₁₈H₂₃NO₃

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to 300°C

(5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN/SIM

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Remarks

Derivatization


Silylation (TMS) with BSTFA/TMCS
(Reagent: Fluka 15238)

100µl of the hydrochloride standard
(Lilly Research Laboratories),
concentration 1.6mg/ml in methanol,


is evaporated with a gentle flow of nitrogen. To the residue, 50µl derivatization reagent and 125µl pyridine are added and the mixture incubated for 30min at 60°C. Gentle evaporation with nitrogen is repeated and the residue redissolved in chloroform. The solution is ready GC/MSD analysis.

Results


Derivatization is recommended due to the polar nature of the analyte. The TMS reaction leads to the tri-derivative. The underivatized analyte shows no protonated molecule or molecular adducts in PCI/CH₄. The PCI/NH₃ reaction of the silylated analyte is the method of choice. In PCI/NH₃ SIM the signal/noise ratio for 100pg of analyte is calculated as approximately 30/1.




EI-Spectrum, Ractopamine, underivatized, molecular mass = 301 u




PCI/CH₄-Spectrum, Ractopamine, underivatized, molecular mass = 301 u



EI-Spectrum, Ractopamine, TMS derivative, m/z 517; M⁺




PCI/CH₄-Spectrum, Ractopamine, TMS derivative, m/z 518, 546, 558; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺



PCI/NH₃-Spectrum, Ractopamine, TMS derivative, m/z 518, 535; [M + H]⁺, [M + NH₄]⁺

PCI/NH₃– SIM Mode



**Ractopamine, TMS derivative, 100pg,
Retention Time: 11.58min
Ions: m/z 428, 518; Signal/Noise \approx 30/1**

Steroids

Cholesterol CAS-Nr. 57-88-5

Molecular Formula: C₂₇H₄₆O

Estradiol CAS-Nr. 7681-76-7

Molecular Formula: C₁₈H₂₄O₂

Estrone CAS-Nr. 443-48-1

Molecular Formula: C₁₈H₂₂O₂

Testosterone CAS-Nr. 58-22-0

Molecular Formula: C₁₉H₂₈O₂

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 0.7ml/min, 30cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

60°C (1min) – 30°C/min to 270°C

(5.5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: ECNI/CH₄ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage Scan: Tune + 400V

EM Voltage SIM: Tune + 800V

Optimization of Tune Parameters

Emission Current: Tune + 150µA

Flow (Setting): 3ml/min (60)

Remarks

Derivatization

The analytes (Cholesterol/Estradiol/Estrone/Testosterone SIGMA C 8667/E 8750/ E 9750/ T 1500) can be derivatized with BSTFA and with Pentafluorophenyl/DMCS as well as other approaches.

Derivatization with TMS

To 200µl of the analytes at a concentration of 1ng/µl each in chloroform, 100µl of the derivatization reagent (BSTFA, Fluka 15238) is added and the mixture is incubated for 30min at 60°C. After evaporation under a gentle flow of nitrogen, the residue

is redissolved in ethyl acetate; the solution is ready for injection.


Derivatization with PFPh

The analytes, concentration 1-2mg/ml each, are dissolved in 100µl acetonitrile and 100µl pyridine is added. The mixture is shaken (Vortex-Minisher) for 2min and to the clear solution 50µl of derivatization reagent (Pentafluorophenyl/DMCS, Fluka 76750) is added. Be unconcerned with the precipitate. The reaction takes 20min at room temperature. After addition of 800µl chloroform the precipitate dissolves and an aliquot is diluted prior to injection.


Caution: PFPh derivatives are unstable even at low temperatures.

Results


The TMS derivatives show weak response in ECNI Scan mode, e.g. for Cholesterol, concentration 18ng/µl, signal/noise ratio ≈ 10/1. Favourable responses are obtained with the PFPh derivatives of Estrone and Estradiol: low concentrations (≈ 200fg/µl) can be detected in SIM mode.




EI-Spectrum, Cholesterol TMS derivative, m/z 458; M⁺




EI-Spectrum, Cholesterol PFPh derivative, m/z 610; M⁺




ECNI/CH₄-Spectrum, Cholesterol TMS derivative, *m/z* 458; M⁺




ECNI/CH₄-Spectrum, Cholesterol, PFPh derivative, *m/z* 610; M⁺




EI-Spectrum, Estradiol TMS derivative, *m/z* 416; M⁺




EI-Spectrum, Estradiol PFPPh derivative, m/z 720; M^+




ECNI/ CH_4 -Spectrum, Estradiol TMS derivative, m/z 416; M^-




ECNI/ CH_4 -Spectrum, Estradiol PFPPh-derivative, m/z 720; M^-




El-Spectrum, Estrone TMS derivative, m/z 342; M^+




El-Spectrum, Estrone PFPh-derivative, m/z 494; M^+




ECNI/CH₄-Spectrum, Estrone TMS-derivative, M^+ : 342 m/z




ECNI/CH₄-Spectrum, Estrone PFP derivative, m/z 494; M⁺




EI-Spectrum, Testosterone, TMS derivative, m/z 360; M⁺



EI-Spectrum, Testosterone, PFP derivative, m/z 512; M⁺




ECNI/CH₄-Spectrum, Testosterone, TMS derivative, m/z 360; M⁺



ECNI/CH₄-Spectrum, Testosterone, PFPh derivative, molecular mass = 512 u

ECNI/CH₄ – SIM



**Estrone (left); Estradiol (right);
SIM Ions: m/z 494; m/z 720; 0.2pg/μl each**

Tetrahydrocannabinol

(d9-THC) CAS-Nr. 1972-08-3

Molecular Formula: C₂₁H₃₀O₂

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to 300°C

(5min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄ – SCAN

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Emission Current: 150µA

Remarks

Derivatization

a) Trifluoroacetylation (TFA) with MBTFA (Reagent: Fluka 65943)

b) Reaction with Pentafluoropropionic Acid Anhydride (PFPA) (Reagent: Fluka 77292)

a) 100µl of the standard (SIGMA T 4764) at a concentration of 100ng/µl is dissolved in ethyl acetate and evaporated with a gentle flow of nitrogen. To the residue 50µl derivatization reagent is added and the solution is incubated for 30min at 80°C. Gentle evaporation with nitrogen is repeated and the residue redissolved in ethyl acetate.

b) Procedure as described above. After the first evaporation of the residue, 80µl of the derivatization


reagent and 20µl of hexafluoroisopropanol (Fluka 52517) are added and the reaction mixture is incubated for 30min at 70°C. Then evaporation, dilution and GC/MSD analysis.

Results


Even the underivatized analyte at a concentration 10ng/µl can be measured without chromatographic discrimination. The described derivatization reactions lead to THC by-products, not considered here. The trifluoroacetylated analyte can be measured in ECNI mode, however the PFPA derivative shows no significant ECNI spectrum. The cited literature give additional information about THC determinations in PCI and ECNI modes.

Literature


“Cannabinoids in Blood: Advantage of Positive Chemical Ionization in Mass Spectroscopy”
Harry Prest,
Agilent Pub. Nr. 5967-6103
“Determining Cannabinoids in Blood Using Electron Capture Negative Chemical Ionization ...”
Harry Prest,
Agilent Pub. Nr. 5967-6331




EI-Spectrum, Tetrahydrocannabinol, underivatized, m/z 314; M⁺




El-Spectrum, Tetrahydrocannabinol, TFA derivative, m/z 410; M^+




El-Spectrum, Tetrahydrocannabinol, PFPA derivative, m/z 460; M^+



PCI/ CH_4 -Spectrum, Tetrahydrocannabinol, PFPA derivative, m/z 461, 489, 501; $[M + H]^+$, $[M + C_2H_5]^+$, $[M + C_3H_5]^+$



PCI/NH₃-Spectrum, Tetrahydrocannabinol, PFPA derivative, m/z 461, 478; [M + H]⁺, [M + NH₄]⁺



ECNI/CH₄-Spectrum, Tetrahydrocannabinol, TFA derivative, m/z 410; M⁺


THC COOH

11-Nor-d9-Tetrahydrocannabinol

-9-Carboxylic Acid

CAS-Nr. 56354-06-4

Molecular Formula: $C_{21}H_{28}O_4$



GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25 μ m

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

80°C (1min) – 30°C/min to 300°C

(3min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄/NH₃ – SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Buffer Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Emission Current: 150 μ A

Remarks

Derivatization

a) Silylation (TMS) with BSTFA/TMCS (Reagent: Fluka 15238)

b) Reaction with Pentafluoropropionic Acid Anhydride (PFPA) (Reagent: Fluka 77292)

a) 100 μ l of the standard (SIGMA N 3142) at a concentration of 100ng/ μ l is dissolved in ethyl acetate and evaporated with a gentle flow of nitrogen. To the residue, 50 μ l of derivatization reagent is added and the solution is incubated for 30min at 80°C.

Gentle evaporation with nitrogen is repeated and the residue redissolved in ethyl acetate.

b) Procedure as described above. After the first evaporation, 80 μ l of the derivatization reagent and 20 μ l of hexafluoroisopropanol (Fluka 52517) are added to the residue and the reaction mixture is incubated for 30min at 70°C. Then evaporation, dilution and GC/MSD analysis.

Results

Derivatization is recommended for this polar analyte. The analyte response is related to the choice of reagent/buffer gases applied. The PFPA derivative is more suitable for ECNI measurements and the highest sensitivity is achieved in ECNI/NH₃ SIM mode.


Literature

“Zum Nachweis von 11-nor-d9-Tetrahydrocannabinolcarbonsäure in menschlichen Haaren mittels GC/MS”


U. Dressler, Dissertation 2000, Med. Fakultät der Universität München

“Determining Cannabinoids in Blood Using ECNCI with HP 5973 CI MSD”


Harry Prest, Agilent Publication Nr. 5967-6331




EI-Spectrum, Tetrahydrocannabinol Carboxylic Acid, TMS derivative, m/z 488; M^+




El-Spectrum, Tetrahydrocannabinol Carboxylic Acid, PFPA derivative, m/z 640; M^+




PCI/ CH_4 -Spectrum, Tetrahydrocannabinol Carboxylic Acid, TMS derivative, m/z 489, 517, 529; $[M + H]^+$, $[M + C_2H_5]^+$, $[M + C_3H_5]^+$




PCI/ NH_3 -Spectrum, Tetrahydrocannabinol Carboxylic Acid, TMS derivative, m/z 489; $[M + H]^+$




ECNI/CH₄-Spectrum, Tetrahydrocannabinol Carboxylic Acid, TMS derivative, *m/z* 489; [M + H]⁺



PCI/CH₄-Spectrum, Tetrahydrocannabinol Carboxylic Acid, PFPA derivative, *m/z* 641, 669, 681; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺




PCI/NH₃-Spectrum, Tetrahydrocannabinol Carboxylic Acid, PFPA derivative, *m/z* 641, 658; [M + H]⁺, [M + NH₄]⁺




ECNI/CH₄-Spectrum, Tetrahydrocannabinol Carboxylic Acid, PFPA derivative, m/z 620, 640; [M - HF]⁻, [M]⁻

SIM Mode, PFPA derivative

0.5pg/μl each, Ions: m/z 620, 640



ECNI/CH₄, Tetrahydrocannabinol Carboxylic Acid



ECNI/NH₃, Tetrahydrocannabinol Carboxylic Acid

Ion Mode	Derivative	Concent.	Acquis. Mode	Approximate Signal/Noise Ratio
EI	TMS	10ng/μl	Scan	650/1
EI	PFPA	10ng/μl	Scan	620/1
PCI/CH ₄	TMS	10ng/μl	Scan	40/1
PCI/NH ₃	TMS	10ng/μl	Scan	80/1
NCI/CH ₄	TMS	10ng/μl	Scan	10/1
PCI/CH ₄	PFPA	10ng/μl	Scan	140/1
PCI/NH ₃	PFPA	10ng/μl	Scan	75/1
NCI/CH ₄	PFPA	10pg/μl	Scan	15/1
NCI/CH ₄	PFPA	0.5pg/μl	SIM	40/1
NCI/NH ₃	PFPA	0.5pg/μl	SIM	60/1

Table: Tetrahydrocannabinol Carboxylic Acid, Responses

trans Zearalenone

CAS-Nr. 17924-92-4

Molecular Formula: C₁₈H₂₂O₅

GC Parameters

Column: HP-5ms

Agilent Part Nr. 19091S-433

30m x 0.25mm x 0.25µm

Carrier Gas: Helium

Flow: 1.2ml/min, 40cm/sec

Mode: Constant Flow

Injection: Pulsed splitless, 250°C

Oven Temp. Program

70°C (1min) – 25°C/min to 300°C

(4min)

MS Parameters

Mode: EI – SCAN

Tune: Atune

Temperatures:

Source 230°C, Quad 150°C

Mode: PCI/CH₄ – SCAN

Reagent Gas: Methane

Flow (Setting): 1ml/min (20)

Tune: PCI-Methane Autotune

Temperatures:

Source 250°C, Quad 106°C

EM Voltage: Tune + 400V

Mode: PCI/NH₃ – SCAN

Reagent Gas: Ammonia

Flow (Setting): 1ml/min (20)

Tune: PCI-Ammonia Tune File

EM Voltage: Tune + 400V

Mode: ECNI/CH₄– SCAN/SIM

Buffer Gas: Methane

Flow (Setting): 2ml/min (40)

Tune: ECNI-Methane Tune File

Temperatures:

Source 150°C, Quad 106°C

EM Voltage: Tune + 400V

Remarks

Derivatization

a) Silylation (TMS) with BSTFA/TMCS (Reagent: Fluka 15238)

b) Reaction with Pentafluoropropionic Acid Anhydride (PFPA) (Reagent: Fluka 77292)

a) 50µl of the standard (SIGMA Z 2125) at a concentration of 10mg/ml in methanol is evaporated with a gentle flow of nitrogen. To the residue, 50µl of derivatization reagent is added and the solution is incubated for 30min at 60°C.

Gentle evaporation with nitrogen is repeated and the residue dissolved in ethyl acetate.

b) Procedure as described above. After the first evaporation to the residue, 40µl of the derivatization reagent and 20µl of hexafluoroisopropanol (Fluka 52517) are added and the reaction mixture is incubated for 30min at 70°C. Then evaporation, dilution and GC/MSD analysis.

Results

Underivatized analyte:

In EI scan mode measurements in the concentrations range of 1ng to 10ng are achievable. In PCI mode, ammonia is the preferred reagent gas.


TMS derivative:

The analyte reacts to make the mono and di-derivative with a response ratio of ≈ 0.4 : 1. PCI/NH₃ responses are a factor of three more sensitive compared to the PCI/CH₄.


PFPA derivative:

The di-derivative is preferentially produced (90%) and elutes prior to the mono-derivative. PCI/CH₄ is the preferred mode compared to ammonia reagent gas.


In ECNI/CH₄ mode the mono-derivative dominates by a factor of more than 10 to one. In SIM a concentration of 0.5pg/µl produces approximately 30:1 signal:noise. In addition to the mentioned derivatization reactions, pentafluorophenyl reaction was tested and resulted in less sensitive response when compared to the PFPA data.




EI-Spectrum, Zearalenone, underivatized, m/z 318; M⁺




PCI/CH₄-Spectrum, Zearalenone, underivatized, *m/z* 319, 347, 359; [M + H]⁺, [M + C₂H₅]⁺, [M + C₃H₅]⁺




PCI/NH₃-Spectrum, Zearalenone, underivatized, *m/z* 319, 336; [M + H]⁺, [M+NH₄]⁺




EI-Spectrum, Zearalenone, TMS mono-derivative, *m/z* 390; M⁺




EI-Spectrum, Zearalenone, TMS di-derivative, m/z 462; M^+




PCI/ CH_4 -Spectrum, Zearalenone, TMS mono-derivative, m/z 391, 419, 431; $[M + H]^+$, $[M + C_2H_5]^+$, $[M + C_3H_5]^+$




PCI/ CH_4 -Spectrum, Zearalenone, TMS di-derivative, m/z 463, 491, 503; $[M + H]^+$, $[M + C_2H_5]^+$, $[M + C_3H_5]^+$




PCI/NH₃-Spectrum, Zearalenone, TMS mono-derivative, m/z 391, 408; $[M + H]^+$, $[M+NH_4]^+$




PCI/NH₃-Spectrum, Zearalenone, TMS di-derivative, m/z 463; $[M + H]^+$




EI-Spectrum, Zearalenone, PFPA mono-derivative, m/z 464; M^+




EI-Spectrum, Zearalenone, PFPA di-derivative, m/z 610; M^+




PCI/ CH_4 -Spectrum, Zearalenone, PFPA mono-derivative, m/z 465, 493; $[M + H]^+$, $[M + C_2H_5]^+$




PCI/ CH_4 -Spectrum, Zearalenone, PFPA di-derivative, m/z 611, 639, 651; $[M + H]^+$, $[M + C_2H_5]^+$, $[M + C_3H_5]^+$




PCI/NH₃-Spectrum, Zearalenone, PFPA mono-derivative, *m/z* 465, 482; [M + H]⁺, [M + NH₄]⁺



PCI/NH₃-Spectrum, Zearalenone, PFPA di-derivative, *m/z* 628; [M + NH₄]⁺




ECNI/CH₄-Spectrum, Zearalenone, PFPA mono-derivative, *m/z* 316, 444, 464; [M - 148(PFPA+H)]⁻, [M - HF]⁻, [M]⁻



ECNI/CH₄-Spectrum, Zearalenone, PFPA di-derivative, molecular mass = 610 u m/z 463 ; [M-147]⁻ (PFP = 147 u)

ECNI/CH₄ – SIM Mode




Zearalenone, PFPA mono-derivative, 200fg,
Retention Time: 10.68min, Ion: m/z 316;
Signal/Noise \approx 30/1

The analysis of pharmacologically relevant compounds is carried out using different analytical techniques. The combination of gas chromatography (GC) with mass spectrometry (MS) is one of the most frequently applied because it offers high separating power, which is advantageous in analyzing complex mixtures, and reliable identification of unknown compounds. In the field of human and veterinarian medicine, drugs of abuse can be unambiguously determined and accurately quantitated.

This brochure emphasizes the standard mass spectral technique of electron impact ionization (EI) and both chemical ionization (CI) techniques; positive chemical ionization (PCI) and electron capture negative ionization (ECNI or NCI). The data presented unmistakably indicates CI is not merely a supplement to EI, but for most of the documented examples, improves analytical results related to compound selectivity and detection sensitivity. To assist and encourage the user in exploring CI, this brochure also presents an elementary understanding of CI-theory and useful practical operating advice for CI on the Agilent MSD.

Data for 43 compounds are presented from a wide variety of drug classes; barbiturates, benzodiazepines, β -agonists, narcotics and steroids and others. All GC and MS method parameters, derivatisation techniques, and the resulting EI/PCI/ECNI mass spectra are documented. The results are briefly discussed for the different operating modes and reagent gases applied. Attention is given to the quantitative data obtained especially in terms of signal-to noise. The intention being to assist the user in choosing successful conditions and modes to analyze samples for particular drugs at trace concentrations.



© Copyright Agilent Technologies Company, 2001.
All rights reserved. Reproduction, adaption, or translation
without prior written permission is prohibited, except as
allowed under the copyright laws.

Printed in Germany
February, 2002

Publication Number 5988-3921EN

www.agilent.com/chem



Agilent Technologies