

PMA-12

Photonic multichannel analyzer



Scientific applications

- UV to visible spectroscopy
- Gas chromatography
- Fluorescence spectroscopy
- Raman scattering
- Luminous efficiency measurement
- Discharge spectrum analysis
- Chemiluminescence analysis
- Combustion analysis
- Liquid chromatography
- Micro spectroscopy

Industrial applications

- Water quality testing
- UV radiation measurements
- Evaluation of light emitting devices and light sources
- Plasma monitoring
- Photobiological safety assessment
- Chromaticity measurements
- Combustion monitoring
- Color filter evaluation
- Film thickness measurements

PRODUCT INTRODUCTION

Photo-detector

Spectrometer

Power supply

Use of an optical fiber input
makes spectral measurements easy.

The PMA-12 is a compact spectral measurement system that combines a spectrometer and optical detector into one unit. Because of the high sensitivity, spectra can easily be obtained in many applications, just by bringing the optical fiber close to the sample without the connection to a special light collection system. Since the spectrometer and photo-detector are manufactured with high machine accuracy, the PMA-12 is stable and can be used with confidence for long periods of time. The wavelength axis and spectral response characteristics are already calibrated, so spectral measurements can be carried out easily and accurately.



▲ C14880-01

C14880-01 Low stray light model

This model realizes low stray light and enables highly accurate spectrum analysis by reviewing the optical layout. By using a sensor with a built-in cooling element, low noise and highly reproducible measurements have been achieved.

C14631-01, -02, -03 High sensitivity superior cost-performance model

The most compact high-performance model in the PMA-12 series. Highly accurate spectrum analysis from basics to applications can be realized at a low price. The wavelength range for measurements is 300 nm to 800 nm for the C14631-01 and 250 nm to 840 nm for the C14631-02 and 300 nm to 1040 nm for the C14631-03.

C10028-01, -02 Near infrared model


These are models using InGaAs linear image sensors which are capable of measuring reflection and absorption spectra in the near infrared with a large dynamic range. The wavelength range for measurements is 900 nm to 1650 nm for the C10028-01 and 1600 nm to 2350 nm for the C10028-02.

C10027-01, -02 Ultra-high sensitivity model

The ultra-sensitive model that can measure a wide wavelength range from ultraviolet to near-infrared with high wavelength resolution by combining with a compact Czerny-Turner spectrometer. The wavelength range for measurements is 200 nm to 950 nm for the C10027-01 and 350 nm to 1100 nm for the C10027-02.


Features

- Spectrometer, photo-detector and power supply in a compact unit
- Real-time measurements (Simultaneous measurement of multiple wavelengths possible)
- Easy measurements with optical fiber
- Spectral response and wavelength calibrated
- Support many applications with the option




● Spectral response (Typ.)

BT-CCD linear image sensor




InGaAs linear image sensor




Measurement modes

- Standard measurements**
 This is the most basic measurement mode.
 Applications: e.g. emission spectra for light sources, fluorescence, plasma and etc.
- Reflective measurements**
 This is the measurement mode for finding spectral reflectance.
 Applications: e.g. reflectance measurements for optical filters, coatings and etc.
- Transmittance and absorption measurements**
 This is the measurement mode for finding spectral transmittance and absorption.
 Applications: e.g. measurements of transmittance and absorption in optical filters, films, solutions and etc.
- Chromaticity measurements (light-source color)**
 This is the measurement mode for finding the light-source color for luminous bodies.
 Applications: e.g. color evaluation in light sources for illumination, LEDs and etc.
- Chromaticity measurements (object color)**
 This is the mode for finding the color of objects that are either reflective or transmit light.
 Applications: e.g. color evaluation of paint, fabric, printed matter and etc.




Display modes


Spectrum display




Display of changes over time




3-D display




Reflectivity display




Transmittance display




Absorbance display (OD)



Color coordinate display




Spatial color coordinate display



APPLICATION EXAMPLES


Light source measurements

Measurement of emission spectra in light sources such as lamps and LEDs




Emission spectrum measurements

Emission spectrum measurements for plasma, electric discharge, ablation and the like



Reflective spectrum measurements


Measurement of spectral reflectance in optical filters, anti-reflective films (AR coatings) and the like



APPLICATION EXAMPLES


Object color measurements

Object color measurement of paint, fabric, printed matter and the like




Absorption spectrum measurements

Spectral transmittance and absorption measurements in optical filters, films, solutions and the like



Microscopic spectral measurements


Spectral distribution measurements under a microscope



APPLICATION EXAMPLES

Emission spectrum measurements


For fluorescent samples such as fluorescent lamps and EL devices




< Configuration >

- Standard PMA-12 configuration (C14880, C14631, C10027, etc.)
- Options
- Excitation light source: laser, xenon lamp, etc.
- Sample Holder for transmission and fluorescence measurement A6751


Fluorescence indicator (Fluorescein) emission spectrum



Chemiluminescence emission spectrum



Emission spectrum of fluorescent materials (Fluorescent lamp)




< Applications >

- Fluorescence spectroscopy
- Monitoring chemical light emissions

Film thickness measurements

Film thickness measurements using spectral reflectance or transmittance




< Configuration >

- Standard PMA-12 configuration (C10027)
- Options
- Halogen lamp L6758-11E
- Optical split fiber VIS to NIR 2 m A10193-02
- Film thickness measurement software U10339-01

< Applications >

- Monitoring thin film growth
- Film thickness management
- Resist film thickness measurements

ITO film interference spectrum




Optical Gauge series

C10178, C10323

We can offer a special machine for film thickness measurements. Please refer to the details in a specific brochure.

Quantum yield measurement system

Measurement of quantum yield, external quantum efficiency, brightness light distribution characteristics




< Configuration >

- Standard PMA-12 configuration (C10027)

< Applications >

- Research of fluorescence materials in physics or chemistry
- Quantum yield measurement of emission materials
- Internal quantum yield measurement of fluorescence materials

[Screen showing emission spectrum]



Absolute PL quantum yield spectrometer C9920-02,-02G,-03,-03G

External quantum efficiency measurement system C9920-12

Light distribution measurement system C9920-11

We can offer a special machine for quantum yield measurements. Please refer to the details in a specific brochure.

OPTIONS



**Sample Holder for transmission and fluorescence measurement
A6751**

This is a dedicated holder with an integrated condensing lens for the use with vials.



**Reflection measurement optics
A9665**

These are optics making it possible to illuminate the sample at 45° to the light source and measure the reflected light.



**Optical split fiber
A10193-01,-02**

It is very useful for reflectance measurement or film thickness measurement. We have two kinds of fiber. One is A10193-01 for from UV to visible light and the other is A10193-02 for from visible to NIR light range.



**C-mount fiber adapter
A6399**

This is an adapter for securing the fiber input optics to the C-mount of a microscope or the like. The A6399 is usable in the UV to NIR.



**C-mount adapter for positioning
A9607**

In addition to the function of the C-mount fiber adapter, the measurement position can be checked. The A9607 is usable in the UV to NIR.



**Objective lens 50 mm F3.5
A4869**

Condensing lens for UV. $f=50$ mm, F3.5 (A6399 or A9607 required)



**Attenuation fiber adapter
A10474-01**


This adaptor is used when the light power is too strong. It can reduce the input light power by using a pinhole. (fading rate approx. 1/20 to 1/500)



**Halogen lamp
L6758-11E**

This is a halogen light source with output wavelengths from 400 nm to 1600 nm for excitation and absorption measurements.

■ L6758-11E emission spectrum




* Light guide connector A10194-01 is needed to connect with 2 split fiber.



**Xe light source High stability 150 W
L6759**

This is a high stability xenon light source with output wavelengths from 250 nm to 1600 nm for excitation and absorption measurements.

■ L6759 emission spectrum



Software library U10472-01

This is the software library which controls the PMA-12 series.

Color measurement library U10473-01

This is the software library which controls the PMA-12 series and calculates the chromaticity.

SPECIFICATIONS


Model	C14631-01	C14631-02	C14631-03	C14880-01	C10027-01	C10027-02	C10028-01	C10028-02
Photo-detector	BT- CCD linear image sensor						InGaAs linear image sensor	
Wavelength (nm)	300 to 800	250 to 840	300 to 1040	200 to 990	200 to 950	350 to 1100	900 to 1650	1600 to 2350
Wavelength resolution (FWHM)*1	≤ 3 nm	≤ 3 nm (Less than 750 nm)	≤ 4 nm	< 2.5 nm	< 2 nm	< 2.5 nm	< 9 nm	
Wavelength accuracy	<±0.3 nm		<±0.5 nm	<±0.3 nm	<±0.75 nm		<±3.2 nm	
Exposure time (Internal trigger Mode)	18 ms to 64 s			19 ms to 64 s			5 ms to 64 s	5 ms to 0.05 s
Number of photosensitive device channels	1024 ch						256 ch	
Pixel size	24 μm × 1392 μm			24 μm × 2928 μm			50 μm × 250 μm	
Device cooling temperature	0 °C			-10 °C	-15 °C		-10 °C	
Read-out noise (electrons) (Max.)	16						18 750	
Dark current (electrons/scan) (Max.)	128 (0 °C : 20 ms)			64 (-10 °C : 20 ms)	32 (-15 °C : 20 ms)		163 000 (-10 °C : 20 ms)	6.47 × 10 ⁸ (-10 °C : 20 ms)
AD resolution	16 bit							
Spectrograph	Concave spherical grating type			Czerny-Turner type				
Spectrograph F number	3			4				
Fiber type	Bundled fiber Φ12 mm SUS tube							
Fiber length	2 m			1.5 m				
Fiber receiving area	Φ1 mm							
External trigger input	TTL level/High impedance							
Interface	USB 2.0 ²			USB 3.0 ²	USB 2.0 ²			
Power supply	AC 100 V to AC 240 V, 50 Hz/60 Hz (Power supply voltage variation ±10 %)							
Power consumption	Approx. 40 VA			Approx. 50 VA	Approx. 70 VA		Approx. 50 VA	
Ambient operating temperature	+10 °C to +35 °C			+10 °C to +30 °C				

*1 Confirmed with mercury and argon atomic beams. *2 1.5 m cable is included as standard.


Dimensional outlines (Unit : mm)

● Main unit


C14631-01,-02,-03 (Approx. 2.6 kg)




C14880-01 (Approx. 5.4 kg)




C10027-01, -02 (Approx. 5.7 kg)



C10028-01, -02 (Approx. 9 kg)



● Fiber input optics (Approx. 100 g)



Basic software for PMA-12 U6039-01

- Measurement functions Monitoring measurement
Data measurement
- Temporal resolution measurement functions ... Temporal fluctuation of spectra
Temporal fluctuation in reflectivity and transmissivity
- Data acquisition condition settings Exposure time settings
Memory integration count assignment
- Calibration/correction Wavelength axis calibration
Sensitivity inconsistency calibration
Dark current correction
- Display functions Spectrum display
Display temporal waveform fluctuations
- Wavelength axis display Wavelength, Wavenumber, Raman shift, energy (eV)
- Brightness axis display Linear, Logarithmic
- Cursor functions Wavelength (wavenumber, etc.) vs. intensity
Peak detection
FWHM measurement
Integrated intensity
- Other functions Smoothing
Differential waveform
Color calculation (XYZ, xy, uv, Lab)

- Windows is registered trademark of Microsoft Corporation in the U.S.A.
 - Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.
 - Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
 - Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.
- Specifications and external appearance are subject to change without notice.
© 2022 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-433-8031, E-mail: export@sys.hpk.co.jp

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509 031 00, Fax: (46)8-509 031 01 E-mail: info@hamamatsu.se

Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6, 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201 Tower B, Jiaxing Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 8F-3, No.158, Section 2, Gongdao 5th Road, East District, Hsinchu, 300, Taiwan R.O.C., Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

Cat. No. SDSS0008E20
JUN/2022 HPK
Created in Japan